


Product Catalog 2017



We believe that
*"Mobility is a basic human right that should be
affordable to everyone, everywhere."*








www.swisslimbs.org

For orders or information email: info@swisslimbs.org






Section « T »: Tibial Prostheses

Item ID	Item Name	Image	Item Description	Order specification	Order unit
T.0001	Adult SACH foot		Solid Ankle Cushioned Hill Foot with toes. Polyuretane and wood. Thread M10. Weight limit 90kg to 100kg.	Colour: Flesh or Brown Side: Left or Right. Size: 20cm - 30cm	single piece
T.0002	Adult SACH foot - Long Stump		Solid Ankle Cushioned Hill Foot with toes. Colour flesh. Polyuretane and wood. Thread M10. Weight limit 100kg.	Side: Left or Right. Size: 21cm - 27cm	single piece
T.0003	Children SACH foot		Solid Ankle Cushioned Hill Foot with toes. Colour flesh. Polyuretane and wood. Thread M8. Weight limit 35kg to 45kg.	Side: Left or Right. Size: 12cm - 19cm	single piece
T.0004	Children SACH foot		Solid Ankle Cushioned Hill Foot. Colour flesh. Polyuretane and wood. Thread M8. Weight limit 35kg.	Side: Left or Right. Size: 12cm - 17cm	single piece
T.0005	Single Axis Foot One-Hole		Heel height: 10mm Weight Limit: 80kg - 90kg Wooden Keel	Colour: Flesh Side: Left or Right. Size: 21 - 27cm	single piece
T.0006	Single Axis Foot Two-Holes		Heel height: 10mm Weight Limit: 80kg - 90kg Wooden Keel	Colour: Flesh Side: Left or Right. Size: 21 - 30cm	single piece
T.0007	Single Axis Foot Adapter One Hole		Adapter with one screw Stainless Steel weight 264-286gr Titanium weight 187-192 gr Weight Limit: 100kg	Stainless steel or titanium Size: 22 - 25cm, or Size: 26 - 27cm	single piece

For orders or information email: info@swisslimbs.org

T.0008	Single Axis Foot Adapter Two Holes		Adapter with two screws Stainless Steel weight 320-330gr Titanium weight 188-197 gr Weight Limit: 100kg	Stainless steel or titanium Size: 22 - 25cm Size: 26 - 30cm	single piece
T.0009	Sach Foot Adaptor		Material: Stainless Steel/Titanium Weight: 40g/60g Weight Limit: 100kg child: 65gr in stainless steel and 60Kg weight limit	Adult or Child Stainless Steel or Titanium	single piece
T.0010	Connection Cap		Connection cap for single axis foot	Sizes from 22-28 cm	single piece
T.0011	Premium High Density Polypropylene Sheets		High Density Co-Polymer. Special formula - optimal no-shrinking properties, highly flexible Sheets Dimensions: 2440x1220mm Monolimb: One sheet of 5mm makes 12 prostheses	Color: Brown/White. Thickness: 2.0 to 3.5 mm 5mm 6mm 9mm 12mm 15mm Note: suggested for Monolimb prostheses are 5 - 6mm	Sheet
T.0012	Adaptor - high resistance		Adaptor for SACH Foot. Child: diam. 4cm; height 8cm Adult small: diam. 6cm; height 9cm Adult medium: diam. 6cm; height 10cm Adult big: diam. 6cm; height 11cm Thread for adult M10 and for child M8.	Child Adult small Adult medium Adult big	single piece
T.0013	Dummy - reusable		Dummy for thermomolding - reusable (approx 10 times)	Child Adult	single piece
T.0014	Screw for SACH Foot		Screw with washer and nut in steel. Child Foot: diameter M8; length 8cm Adult Foot: diameter M10; length 9cm	M8 M10	Box with 100 pieces per each size and each diameter

For orders or information email: info@swisslimbs.org








T.0015	Soft Liner		EVA Foam with silvershield (antimicrobial lining). Dense foam with smooth surface finish. For all applications in prosthetics and orthotics, including insoles. Sheet size 1050x1050mm	Colour: Flesh or Terra Thickness: 2mm - 35mm <small>Note: suggested size for Trans-Tibial soft liner prostheses are 5 - 7mm suggested size for cosmetic prostheses covering are 3mm</small>	Sheet
T.0016	EVA - Keasy		EVA Foam For prosthetics and orthotics soft sockets.	S: 13X45cm M: 15X45cm L: 17X45cm	single piece
T.0017	Pedilen Rigid for Foam 200		For fabricating compensations and fillings. Fast and safe processing with hand mixing procedure. No Shrinkage. Easy to refinish for sanding and lamination, Penetration Hardness approx 12.50 N/mm2 Drums of 4Kg		single piece
T.0018	Hardener for Pedilen		For Pedilen Duplicating Foam to achieve Pedilen foam in exact corresponding measure. Drums of 4Kg		single piece
T.0019	Plastazote		Plastazote foam: polyethylene foam. Evazote polyethylene/eva copolymer foam. Good chemical resistance. Will thermoform into simple or complex shapes. Applications: Light Splinting, Collars, Padding, Wristsplints, Shoe build up's and Lower limb prosthetic cosmetic finishing Sheet size 1000x750mm or 1500mmx1000mm	Thickness: 2mm - 27mm sheet small: 1000x750mm sheet medium: 1500mmx1000mm colour: white / blue / black / green / red / yellow / flesh <small>Note: suggested thickness for Trans Tibial stump padding is 12mm</small>	Sheet

*Prices will be quoted in CHF converted in USD at current rate and include FOB loading









*Prices do not include transportation costs, insurance, taxes and duties

For orders or information email: info@swisslimbs.org












Section « M »: Ortho-Prosthetic Manufacturing

Item ID	Item Name	Image	Item Description	Order specification	Order unit
M.0001	Rasp Head Fir Cone		Rasp head " Fir Cone " for Inclined arm machine for processing plastics and other hard materials. Fine rasp stroke. Thread for M16 and 5/8"	Note: More rasp heads and polishing tools for inclined arm milling machines are available on request.	single piece
M.0002	Rasp Head Pine Cone		Rasp head " Pine Cone " for Inclined arm machine for processing wood, pedilene, foam, etc. Medium rasp stroke. Thread for M16 and 5/8"	Note: More rasp heads and polishing tools for inclined arm milling machines are available on request.	single piece
M.0003	Rasp Head Spiral-Toothed		Rasp head " Spiral-toothed " for Inclined arm machine for processing plastics and light materials. Medium rasp stroke. Thread for M16 and 5/8"	Note: More rasp heads and polishing tools for inclined arm milling machines are available on request.	single piece
M.0004	Polishing Wheel		Untreated cotton cloth folded in waves, for working plastics and smoothening surfaces. For Inclined arm machine to be mounted on connector.		single piece
M.0005	Sanding Sleeve long		Sanding Sleeve for inclined arm machine to be mounted on sanding drum. Width 200 mm and diameter 72.5mm	Grain small: 50 Grain medum: 80	single piece
M.0006	Sanding Sleeve short		Sanding Sleeve for inclined arm machine to be mounted on sanding holder. Width 70mm and diameter 45mm	Grain small: 40 Grain medum: 80 Grain big: 120	single piece
M.0007	Sanding Heads		Sanding head for inclined arm machine to be mounted on head of router machine. Diameter 45mm	Grain small: 40 Grain medum: 80	single piece





For orders or information email: info@swisslimbs.org

M.0008	Transfer Paper Blue Camouflage		For Design and Colour on all Orthosis and Porstheses Applications: On all thermoplastic materials and also used on polyethylene foams		per sqm
M.0009	Transfer Paper Brown Camouflage		For Design and Colour on all Orthosis and Porstheses Applications: On all thermoplastic materials and also used on polyethylene foams		per sqm
M.0010	Transfer Paper Football BW		For Design and Colour on all Orthosis and Porstheses Applications: On all thermoplastic materials and also used on polyethylene foams		per sqm
M.0011	Transfer Paper Dinosaurs		For Design and Colour on all Orthosis and Porstheses Applications: On all thermoplastic materials and also used on polyethylene foams		per sqm
M.0012	Transfer Paper Blue Jeans		For Design and Colour on all Orthosis and Porstheses Applications: On all thermoplastic materials and also used on polyethylene foams		per sqm
M.0013	Transfer Paper Flowers		For Design and Colour on all Orthosis and Porstheses Applications: On all thermoplastic materials and also used on polyethylene foams		per sqm
M.0014	Transfer Paper Spring Butterflies		For Design and Colour on all Orthosis and Porstheses Applications: On all thermoplastic materials and also used on polyethylene foams		per sqm
M.0015	Transfer Paper Skulls		For Design and Colour on all Orthosis and Porstheses Applications: On all thermoplastic materials and also used on polyethylene foams		per sqm

For orders or information email: info@swisslimbs.org

M.0016	Files Holders for Plaster Flat		Surform Rasp Holder (without blade) Length 250mm		single piece
M.0017	Files Holders for Plaster Round		Surform Rasp Holder (without blade) Length 250mm		single piece
M.0018	Blades for Plaster		Surform blades for working plaster: Length: 250mm	Flat or half round or round Grain: fine or medium Blades: flat, half round, round	single piece
M.0019	Heat Protective Gloves		Isothermal. Special Textile gloves for high heat work (for thermomolding)		single piece
M.0020	Protective Masks		Filtering mask for protection of dusty work in all industry. With minimal breathing resistance and maximal filter performance		single piece
M.0021	Protective Goggles		Plastic Glasses, colourless. Nylon, with folding side protector. Crystal clear feature		single piece
M.0022	Plaster Mixing Bowl		Rubber black Capacity: 0,5Kg Diameter: 7/12cm		single piece
M.0023	Plaster Smoothing Tool		Spatula for plaster work Blade width: 12mm		single piece
M.0024	Glue Dispenser		Glue Container for adhesive Capacity: 0.45L Empty Comes with one brush		single piece
M.0025	Pensils for Plaster and Plastic		blue for plaster red for plastic and metal		single piece
M.0026	Smoothing/ Buffing polisher		Finishing and polishing of varnished surfaces of prosthesis sockets and all plastics		single piece

For orders or information email: info@swisslimbs.org





M.0027	Caliper		to 1-20mm, capacity 300mm inside measure and long jaws 90mm		single piece
M.0028	Bubble Level		Aluminium level, vertical and horizontal measurement		single piece
M.0029	Taylor's Scissors		Length 250mm, special trimming scissor		single piece
M.0030	Tool Box		with over 100 tools for professional use		single piece

*Prices will be quoted in CHF converted in USD at current rate and include FOB loading

*Prices do not include transportation costs, insurance, taxes and duties

For orders or information email: info@swisslimbs.org



Section « X »: Ortho-Proshtetic Machines

Item ID	Item	Image	Item Description	Order unit
X.0001	Bench Finishing Machine		<p>Two-phase speed for the model 1st speed 1400g / m - 2nd speed 2800 g / m Engine vertical adjustable speed 2800g / m Silent fan internal suction of fine dust Dust recovery system with removable trays Internal filter fabric and back closure with filtering cells Single-phase 230V - 50Hz TRI-phase: 400V - 50Hz</p> <p>Dimensions and weight Width: 8000 mm - Depth: 760 mm Height: 1650 mm - Weight: 265 kg Weight with packaging: 290 kg</p>	single piece
X.0002	Inclined Arm Router		<p>Horizontal engine with abrasive pad 75 x 200 mm. M16 thread and cone abrasive T22 in masthead. Two speeds of the crankshaft • 1 speed 1400 rpm/min - 2nd speed 2800 g/min Drive with foot switch. Connection for centralized suction system with tube diameter 150mm on the side. Suction channel adjacent the abrasive pad. Suction mouth adjustable. Single phase: 230 v-50 hz - Tri phase: 400v – 50 hz</p> <p>Dimensions and weight Width: 1100 mm - Depth: 600 mm Height: 1500 mm - Length of motor arm: 500 mm Weight with packing: 170 Kg (gross)</p>	single piece
X.0003	Air Treatment System (For Router)		<p>Powerful Air Treatment for Inclined Arm Router with long-lasting filter in cellulose M CEL with large (12msq) surface Single phase: 230V - 50Hz Aspiration Power: 1450 cubic meters per hour</p> <p>Dimensions and weight Width: 560 mm - Depth: 380mm - Height: 480mm Weight: 115 Kg (with packing: 125 Kg)</p>	single piece
X.0004	Portable Finishing Machine		<p>Small portable compact cutter suction ring and abrasive. Useful for minor adjustments. • Tree full engine bayonet. • Supplied abrasive ring 100x100mm diameter. • Integrated on the motor inlet. • External filter fabric. • Optional metal base with shelves and drawers</p> <p>Single phase: 230V - 50Hz</p> <p>Dimensions and weight Width: 540 mm - Depth: 410 mm Height: 440 mm Weight: 24 kg - Weight with packing: 30 Kg</p>	single piece

For orders or information email: info@swisslimbs.org

X.0005	Flat Oven		<p>Suitable for all types of thermoplastic materials used in orthopedics (polyethylene, polypropylene, PETG etc.) Digital electronic thermostat. Removable base Upper plate self-regulating in aluminum alloy to keep uniform the thickness of the thermoplastic sheets. Dimensions: mm.1250 x 1050 x 860h Usable area 1100 x 865 mm Weight Kg: 135 Heating elements: 2 x 2000 W 220V single-phase or 380V tri-phase</p>	single piece
X.0006	Vacuum Unit		<p>6-Way Vacuum pump for thermoforming. Tub auto model membrane. lever control vacuum. internal pump with automatic shut off. Single phase: 230V - 50Hz Dimensions and weight Width: 500 mm -Depth: 400 mm Height: 500 mm - Weight: 40 kg Weight with packing: 40 kg</p>	single piece
X.0007	Air Compressor		<p>Air Compressor 40L Piston Compressor, fully automatic, via pressure regulator. Intake filter with noise damper. Tank water drain valve. On wheels. Dimensions and weight: 410mm x 970mm x 770mm Max pressure 10 bar</p>	single piece
X.0008	Electric Jig Saw		<p>Jig Saw with 4-stage pendulum stroke and knob handle. Quick clamping system for changing saw blades with no need for tools. Includes: VarioBoX, saw blade, chip guard, sole plate. Speed: 450-3,00 min-1. Weight: 2,3 Kg</p>	single piece
X.0009	Hot-air Gun		<p>Hot-air Gun with high heating power. With stable blower, adjustable temperature, with electronic temperature control. Suitable for wide range of applications. Includes: Adaptors. Temperature Range: 90-600 degrees Celsius. Weight: 0,6 Kg</p>	single piece
X.0010	Cast Cutter		<p>Oscillating blade for cutting plaster. With onroller for continuously variable pre-selection of the oscillating frequency. Heavy duty, robust design. Includes: saw blades diameter: 50 and 65 mm</p>	single piece
X.0011	Hand Drill		<p>Hand Drill with anti-vibration system for easy and fatigue-free use. Electronic pre-selection of speed: Equipment included: transport case, hand grip with depth stop, quickaction chuck. Drilling performance in concrete/steel/wood (22/13/40). Speed: 0-3,200 min-1. Weight: 2,9 Kg</p>	single piece

For orders or information email: info@swisslimbs.org













X.0012	Laser Aligner		Niveau laser for the perfect alignment. Vertical or Horizontal or both options uses 4 common AA Batteries Includes: trypod PLL360	single piece
X.0013	Casting Frame for Above Knee		For functionbal casting under weight	single piece

*Prices will be quoted in CHF converted in USD at current rate and include FOB loading








*Prices do not include transportation costs, insurance, taxes and duties

For orders or information email: info@swisslimbs.org

Section « OR »: **Ready-to-Wear Orthoses**

Item ID	Item	Image	Order specification	Order unit												
OR.0001	Thumb Support		<table><tr><th>Size</th><th>Wrist circumference</th></tr><tr><td>S</td><td>13-15 cm</td></tr><tr><td>M</td><td>15-17 cm</td></tr><tr><td>L</td><td>17-19 cm</td></tr><tr><td>XL</td><td>19-21 cm</td></tr></table> <p>Available LEFT or RIGHT</p> 	Size	Wrist circumference	S	13-15 cm	M	15-17 cm	L	17-19 cm	XL	19-21 cm	single piece		
Size	Wrist circumference															
S	13-15 cm															
M	15-17 cm															
L	17-19 cm															
XL	19-21 cm															
OR.0002	Wrist + Thumb Support		<table><tr><th>Size</th><th>Wrist circumference</th></tr><tr><td>S</td><td>13-15 cm</td></tr><tr><td>M</td><td>15-17 cm</td></tr><tr><td>L</td><td>17-19 cm</td></tr><tr><td>XL</td><td>19-21 cm</td></tr></table> <p>Available LEFT or RIGHT</p> 	Size	Wrist circumference	S	13-15 cm	M	15-17 cm	L	17-19 cm	XL	19-21 cm	single piece		
Size	Wrist circumference															
S	13-15 cm															
M	15-17 cm															
L	17-19 cm															
XL	19-21 cm															
OR.0003	Wrist Splint		<table><tr><th>Size</th><th>Wrist circumference</th></tr><tr><td>XS</td><td>11-13 cm</td></tr><tr><td>S</td><td>13-15 cm</td></tr><tr><td>M</td><td>15-17 cm</td></tr><tr><td>L</td><td>17-19 cm</td></tr></table> <p>Available LEFT or RIGHT</p> 	Size	Wrist circumference	XS	11-13 cm	S	13-15 cm	M	15-17 cm	L	17-19 cm	single piece		
Size	Wrist circumference															
XS	11-13 cm															
S	13-15 cm															
M	15-17 cm															
L	17-19 cm															
OR.0004	Wrist Hand Splint		<table><tr><th>Size</th><th>Wrist circumference</th></tr><tr><td>XS</td><td>11-13 cm</td></tr><tr><td>S</td><td>13-15 cm</td></tr><tr><td>M</td><td>15-17 cm</td></tr><tr><td>L</td><td>17-19 cm</td></tr></table> <p>Available LEFT or RIGHT</p> 	Size	Wrist circumference	XS	11-13 cm	S	13-15 cm	M	15-17 cm	L	17-19 cm	single piece		
Size	Wrist circumference															
XS	11-13 cm															
S	13-15 cm															
M	15-17 cm															
L	17-19 cm															
OR.0005	Knee Support Padded		<table><tr><th>Size</th><th>Knee circumference</th></tr><tr><td>S</td><td>30-34 cm</td></tr><tr><td>M</td><td>34-38 cm</td></tr><tr><td>L</td><td>38-42 cm</td></tr><tr><td>XL</td><td>42-46 cm</td></tr><tr><td>XXL</td><td>46-50 cm</td></tr></table> 	Size	Knee circumference	S	30-34 cm	M	34-38 cm	L	38-42 cm	XL	42-46 cm	XXL	46-50 cm	single piece
Size	Knee circumference															
S	30-34 cm															
M	34-38 cm															
L	38-42 cm															
XL	42-46 cm															
XXL	46-50 cm															
OR.0006	Knee Ranger		<table><tr><th>Size</th><th>Knee circumference</th></tr><tr><td>S</td><td>30-34 cm</td></tr><tr><td>M</td><td>34-38 cm</td></tr><tr><td>L</td><td>38-42 cm</td></tr><tr><td>XL</td><td>42-46 cm</td></tr><tr><td>XXL</td><td>46-50 cm</td></tr></table> 	Size	Knee circumference	S	30-34 cm	M	34-38 cm	L	38-42 cm	XL	42-46 cm	XXL	46-50 cm	single piece
Size	Knee circumference															
S	30-34 cm															
M	34-38 cm															
L	38-42 cm															
XL	42-46 cm															
XXL	46-50 cm															

For orders or information email: info@swisslimbs.org

OR.0007	Hinged Knee Brace		<table><tr><th>Size</th><th>Knee circumference</th></tr><tr><td>S</td><td>30-34 cm</td></tr><tr><td>M</td><td>34-38 cm</td></tr><tr><td>L</td><td>38-42 cm</td></tr><tr><td>XL</td><td>42-46 cm</td></tr><tr><td>XXL</td><td>46-50 cm</td></tr></table> <p>How to measure</p> 	Size	Knee circumference	S	30-34 cm	M	34-38 cm	L	38-42 cm	XL	42-46 cm	XXL	46-50 cm	single piece		
Size	Knee circumference																	
S	30-34 cm																	
M	34-38 cm																	
L	38-42 cm																	
XL	42-46 cm																	
XXL	46-50 cm																	
OR.0008	Ultra Guard Hip Orthosis		<p>Individual size</p> <p>Available LEFT or RIGHT</p>	single piece														
OR.0009	Knee Immobilizer No-Hinged		<table><tr><th>Size</th><th>Length of splint</th></tr><tr><td>S</td><td>40 cm</td></tr><tr><td>M</td><td>50 cm</td></tr><tr><td>L</td><td>60 cm</td></tr><tr><td>XL</td><td>70 cm</td></tr><tr><td>XXL</td><td>70 cm</td></tr></table> <p>How to measure</p> 	Size	Length of splint	S	40 cm	M	50 cm	L	60 cm	XL	70 cm	XXL	70 cm	single piece		
Size	Length of splint																	
S	40 cm																	
M	50 cm																	
L	60 cm																	
XL	70 cm																	
XXL	70 cm																	
OR.0010	Lumber Support		<table><tr><th>Size</th><th>Waist circumference</th></tr><tr><td>S</td><td>65-75 cm</td></tr><tr><td>M</td><td>75-85 cm</td></tr><tr><td>L</td><td>85-97 cm</td></tr><tr><td>XL</td><td>97-110 cm</td></tr><tr><td>XXL</td><td>110-125 cm</td></tr><tr><td>XXXL</td><td>125-145 cm</td></tr></table> <p>How to measure</p> 	Size	Waist circumference	S	65-75 cm	M	75-85 cm	L	85-97 cm	XL	97-110 cm	XXL	110-125 cm	XXXL	125-145 cm	single piece
Size	Waist circumference																	
S	65-75 cm																	
M	75-85 cm																	
L	85-97 cm																	
XL	97-110 cm																	
XXL	110-125 cm																	
XXXL	125-145 cm																	

*Prices will be quoted in CHF converted in USD at current rate and include FOB loading








*Prices do not include transportation costs, insurance, taxes and duties

For orders or information email: info@swisslimbs.org

Section « C »: **Modular & AK Components**

Item ID	Item Name	Image	Description - Specifications	Order specification	Order unit
C.0001	4-Ear Rotable Lamination Adapter - receiver		A) Material: Titanium/Stainless Steel B) Weight: 130g / 90g C) Weight Limit: 100kg Child: A)Steel B) 68g C) 60Kg	Adult or Child Stainless Steel or Titanium	single piece
C.0002	4-Ear Rotable Adapter - receiver		A) Material: Stainless Steel B) Weight: 100g / 90g C) Weight Limit: 100kg Child: A)Steel B) 50g C) 60Kg	Adult or Child	single piece
C.0003	4-Ear Rotable Lamination Adapter		Material: Stainless Steel/Titanium Weight: 200g/150g Weight Limit: 100kg	Length: 45mm/50mm/70mm Stainless Steel or Titanium	single piece
C.0004	4-Ear Lamination Adapter		Material: Stainless Steel/Titanium Weight: 85g/50g Weight Limit: 100kg child: 50gr in stainless steel and 50Kg weight limit	adult or child with or without hole Stainless Steel or Titanium	single piece
C.0005	4-Ear Rotable anchor with pyramid adapter		Material: Stainless Steel/Titanium Weight: 125g/150g Weight Limit: 100kg	Stainless Steel or Titanium	single piece
C.0006	4-Hole Pyramid Receiver Adapter		Material: Stainless Steel/Titanium Weight: 90g/50g Weight Limit: 100kg	Stainless Steel or Titanium	single piece
C.0007	4-Hole Pyramid Adapter		Material: Stainless Steel/Titanium Weight: 100g/50g Weight Limit: 100kg	with or without hole Stainless Steel or Titanium	single piece








For orders or information email: info@swisslimbs.org

C.0008	4-Hole Rotatable Pyramid Adapter		Material: Stainless Steel Weight: 170g Weight Limit: 100kg		single piece
C.0009	4-Hole Rotatable Pyramid Receiver Adapter, adj. Base		Material: Stainless Steel Weight: 180g Weight Limit: 100kg		single piece
C.0010	4-Hole Rotatable Pyramid Receiver Adapter		Material: Stainless Steel/Titanium Weight: 120g/85g Weight Limit: 100kg	Stainless Steel or Titanium	single piece
C.0011	4-Hole Pyramid Receiver Round		Material: Stainless Steel/Titanium Weight: 130g/85g Weight Limit: 100kg	Stainless Steel or Titanium	single piece
C.0012	4-Hole Rotatable Pyramid Receiver Cross		Material: Stainless Steel/Titanium Weight: 100g/80g Weight Limit: 100kg	Stainless Steel or Titanium	single piece
C.0013	3-Prong Pyramid Receiver, Rotatable		Material: Stainless Steel/Titanium Weight: 150g/130g Weight Limit: 100kg	Stainless Steel or Titanium	single piece
C.0014	3-Prong Rotatable Pyramid Adapter		Material: Stainless Steel/Titanium Weight: 160g/130g Weight Limit: 100kg	Stainless Steel or Titanium	single piece

For orders or information email: info@swisslimbs.org

C.0015	3-Prong Pyramid Adapter		Material: Stainless Steel Weight: 120g Weight Limit: 100kg	with or without hole	single piece
C.0016	Tube Clamp Adapter		Material: Stainless Steel/Titanium/aluminium Weight: 85g/70g/80g Weight Limit: 100kg child: 65gr in stainless steel and 60Kg weight limit	Adult or Child Stainless Steel or Titanium or Aluminium	single piece
C.0017	Tube Clamp Adapter Pyramid		Material: Stainless Steel/Titanium Weight: 110g/70g Weight Limit: 100kg child: 55gr in stainless steel and 60Kg weight limit	Adult or Child Stainless Steel or Titanium	single piece
C.0018	Lenght Adjustment Adapter		Material: Stainless Steel/Titanium Weight: 240g/140g Length: 60-75 mm Weight Limit: 85kg	Stainless Steel or Titanium	single piece
C.0019	Double Adapters with Pyramid		Material: Stainless Steel/Titanium Weight: 110-170g/70-115g Weight Limit: 100kg child: 55-70gr in stainless steel and 60Kg weight limit	Adult or Child Stainless Steel or Titanium Length: 32mm Length: 45mm Length: 60mm Length: 75mm	single piece
C.0020	Short Tube Adapter		Material: Stainless Steel/Titanium Weight: 100-120g/60-80g Weight Limit: 100kg child: 55-70gr in stainless steel and 60Kg weight limit	Adult or Child Stainless Steel or Titanium Length: 45mm Length: 60mm Length: 75mm	single piece
C.0021	Double Adapters		Material: Stainless Steel/Titanium Weight: 140-180g/90-120g Weight Limit: 100kg	Length: 32mm Length: 45mm Length: 60mm Length: 75mm	single piece

For orders or information email: info@swisslimbs.org

C.0022	Tube Adapter		Material: Stainless Steel/Titanium/aluminium Weight: 220-310g/150-240g/160-250g Weight Limit: 100kg Diameter: 30mm child: 110-170gr in stainless steel and 60Kg weight limit	Adult or Child Lenght: 220mm Lenght: 420mm	single piece
C.0023	Tube		Material: Aluminium/Carbon Fiber Weight: 130-250g/90-160g Weight Limit: 100kg child: 60-110gr in aluminium and 60Kg weight limit, 145 or 310 mm Length	Adult or Child Aluminium or Carbon Fiber Lenght: 200mm Lenght: 400mm	single piece
C.0024	PP Plate Rotable Adapter		Material: Stainless Steel/Titanium Weight: 100-120g/60-80g Weight Limit: 100kg	Stainless Steel or Titanium	single piece
C.0025	Adapter for Modular BK Thermoplastic Socket - Receiver		Material: Stainless Steel/Titanium Weight: 200g/130g Weight Limit: 100kg child: 140gr in stainless steel and 60Kg weight limit	Adult or Child Stainless Steel or Titanium	single piece
C.0026	Adapter for Modular BK Thermoplastic Socket - Pyramid		Material: Stainless Steel/Titanium Weight: 140g/90g Weight Limit: 100kg	Stainless Steel or Titanium	single piece
C.0027	Sliding Adaptor, Double Clamp		Material: Stainless Steel/Titanium Weight: 250-320g/150-190g Weight Limit: 85kg	Stainless Steel or Titanium Length: 45mm Length: 50mm Length: 60mm Length: 70mm	single piece
C.0028	Sliding Adaptor, Single Clamp		Material: Stainless Steel/Titanium Weight: 300g/180g Weight Limit: 85kg Length: 65mm	Stainless Steel or Titanium	single piece

For orders or information email: info@swisslimbs.org







C.0029	Child Sliding Adaptor, Double Clamp		Material: Stainless Steel/Titanium Weight: 200g/120g Length: 65-75mm Weight Limit: 60kg	Stainless Steel or Titanium	single piece
C.0030	Child Sliding Adaptor, Single Clamp		Material: Stainless Steel/Titanium Weight: 200g/120g Length: 80-90mm Weight Limit: 60kg	Stainless Steel or Titanium	single piece
C.0031	Shuttle Lock, with Pyramide		Weight Limit: 100kg	single hole or double hole on ears	single piece
C.0032	Shuttle Lock, with Clamp Adapter		Weight Limit: 100kg	single hole or double hole on ears	single piece

*Prices will be quoted in CHF converted in USD at current rate and include FOB loading







*Prices do not include transportation costs, insurance, taxes and duties

For orders or information email: info@swisslimbs.org

Section « K »: **Trans Femoral Knee & Hip Joints**

Item ID	Item Name	Image	Description - Specifications	Order specification	Order unit
K.0001	Single Axis Knee Joint With Weight-Activated Brake		"3R15" Single Axis with Weight-Activated Brake Knee flexion angle approx. 150° Height: 40mm (knee center to upper edge) Material: Stainless steel Weight: 520g Weigh Limit: 100kg With Extension Assist	With Extension assist or without extension assist	single piece
K.0002	Single Axis Knee Joint with Manual Lock (Light Weight)		"3R40" Single Axis with Manual Lock Knee flexion angle approx. 140° Height: 40mm (knee center to upper edge) Material: Aluminum Weight: 301g Weigh Limit: 85kg		single piece
K.0003	4-Bar Knee Joint		"3R20" 4-Bar Knee Joint Knee flexion angle approx. 120° Height: 37mm (knee center to upper edge) Material: Stainless steel Weight: 690g Weigh Limit: 100kg		single piece
K.0004	4-Bar Knee Joint		"3R35" 4-Bar Knee Joint Knee flexion angle approx. 120° Height: 30mm (knee center to upper edge) Material: Stainless steel Weight: 690g Weigh Limit: 100kg		single piece
K.0005	Knee Joint for Disarticulation		"3R21" Knee Joint for Disarticulation (for through-knee) Knee flexion angle approx. 110° Height: 18mm (knee center to upper edge) Material: Stainless steel Weight: 900g Weigh Limit: 100kg		single piece
K.0006	Polycentric 4-Bar Knee Joint		"2SR320" Polycentric 4-bar knee Available with built-in proximal pyramid Knee flexion angle approx. 165° Material: Aluminum Alloy Weight: 520g Weigh Limit: 125kg		single piece

For orders or information email: info@swisslimbs.org







K.0007	Pneumatic Knee Joint		"2SR420 Freda" Pneumatic Knee Joint Polycentric 4-bar knee Adjustable Flexion and extension controls Knee flexion angle approx. 145° Material: Aluminum Alloy Weight: 680g Weight Limit: 125kg		single piece
K.0008	All-Terrain Knee		Water Resistant, designed for harsh and wet environments. Automatic mechanical lock, provides excellent stability and smooth transition into swing-phase. Swing-phase control provided by advanced variable friction mechanism and extension assist spring.		Single piece
K.0009	Single Axis Knee Joint		Knee Only ROBUST, Waterproof single axis knee joint		Single piece
K.0010	Child Knee Joint		Child Knee joint with pyramid Weight: 390gr Child Knee joint for disarticulation Weight: 440gr	With or without manual lock with pyramid or for knee disarticulation	single piece
K.0011	4 Bar Hip Joint		For hip disarticulations Stainless Steel weight: 850g Weight limit: 100Kg		single piece
K.0012	Single Axis Hip Joint		For hip disarticulations Stainless Steel weight: 650g Weight limit: 75Kg		single piece

*Prices will be quoted in CHF converted in USD at current rate and include FOB loading

*Prices do not include transportation costs, insurance, taxes and duties

For orders or information email: info@swisslimbs.org

Section « AX »: **Accessories**

Item ID	Item Name	Image	Description - Specifications	Order unit
AX.0001	Valve Seat with Tube		Tube Diameter: Outer: 40mm Inner hole: 24mm	single piece
AX.0002	Valve Seat Ring		Ring Diameter: Outer: 40mm Inner hole: 24mm	single piece
AX.0003	Flat Rubber Valve without Seat		Ring Diameter: Outer: 40mm Inner hole: 24mm	single piece
AX.0004	Automatic Rubber Valve without Seat		Ring Diameter: Outer: 40mm Inner hole: 24mm	single piece
AX.0005	Automatic Valve with Seat Ring		Ring Diameter: Outer: 44mm Inner hole: 31mm	single piece
AX.0006	Stump Sock Cotton		Stump Sock in Cotton Measures: 30cm, 40cm, 50cm	single piece
AX.0007	Stump Sock Nylon		Stump Sock in Nylon Measures: 30cm, 50cm	single piece

For orders or information email: info@swisslimbs.org






AX.0008	PVA Film		PVA Film Thickness: 0.08mm Lenght: 100Mt, Width: 110cm Lenght: 50Mt, Width: 100cm	single piece
AX.0009	Carbon Fiber Stockinette		Carbon Fiber Width: 8cm	single piece
AX.0010	Fiber Glass Stockinette		Fiber Glass Stockinette Colour: White Width: small, medium, large	single piece roll
AX.0011	Cotton Stockinette		Cotton Glass Stockinette Colour: White Width: small, medium, large	single piece roll
AX.0012	Perlon Stockinette		Perlon Stockinette Colour: White Width: small, medium, large	single piece roll
AX.0013	Lamination Bags		Width: 13 to 36cm Length: 50 to 100cm	10 pieces per pack

*Prices will be quoted in CHF converted in USD at current rate and include FOB loading

*Prices do not include transportation costs, insurance, taxes and duties

For orders or information email: info@swisslimbs.org

Section « OJ »: Orthotic Joints


Item ID	Item Name	Image	Description - Specifications	Order unit
OJ.0001	Swiss Lock		Spring Reset Lock Grip Material: Stainless Steel Length: Upper 410mm, Lower 440mm Width: 16mm and 19mm Thickness: 4mm	Pair
OJ.0002	Drop Lock		Dropping Lock Strip Material: Stainless Steel Length: Upper 410mm, Lower 440mm Width: 16mm and 19mm Thickness: 4mm	Pair
OJ.0003	Bending Iron	 	Bending Iron for Orthotics With Flat Tines: mm 4&6, 6&8, 7&9 With Rounded Tines: mm 4&6, 6&8	Single Piece
OJ.0004	Rivet Stick			Single Piece

*Prices will be quoted in CHF converted in USD at current rate and include FOB loading

*Prices do not include transportation costs, insurance, taxes and duties

For orders or information email: info@swisslimbs.org

Section « Z »: **Others**

Item ID	Item Name	Image	Description - Specifications	Order unit
Z.0001	Natural Deer Tallow Balm		100% Natural Deer Tallow Balm 100g - for external use only Protective skin treatment that nourishes skin, heals blisters, calluses and excoriations, protects from irritations	single piece

*Prices will be quoted in CHF converted in USD at current rate and include FOB loading

*Prices do not include transportation costs, insurance, taxes and duties